

PROMOTION

NOV - DEC 2024



PRIME LINE



AeroAcoustic blade



Glass rod fiber optics



Reverse flow protection



Automatic pressure regulation¹



Refined ceramic bearings



Scratch resistance
Iconic Prime body

**25
W**

**Powerful & Precise
Power Head - 25W**

**55
dB**

**Incredibly Quiet
Only 55 dB noise**



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



**BUY 1
GET
25% OFF**
\$1,262.25

**BUY 2
GET
30% OFF**
\$1,178.10
each

**BUY 3
GET
35% OFF**
\$1,093.95
each

Coupling Connection

Kavo
NSK
W&H
Sirona

Power Head 25 Watt

HP21KL
HP21NL
HP21WL
HP21SL

Small Head 18 Watt

HP22KL
HP22NL
HP22WL
- N/A -

\$1,683.00



**Now with 2 Years
FairCare Warranty**

Unconditional warranty
except physical damage

COUPLINGS



✓ MK-dent & Kavo
with LED & water
• **QC6016KW**
\$ 495.00



✓ NSK connection
with LED
• **QC6016NT**
\$ 385.00



✓ W&H connection
with fibre optic
• **QC5016W**
\$ 385.00



✓ Sirona connection
with LED & water
• **QC6016SW**
\$ 495.00

✓ Non-optic
• **QC4014K**
\$ 296.00

✓ Non-optic
• **QC4014NT**
\$ 296.00

✓ Non-optic
• **QC4014W**
\$ 296.00

✓ Non-optic
-NA-

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

ECO LINE TURBINES



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



BUY 2
\$890
each

2 x EcoLine F/Optic



BUY 3+
\$792
each

3 x EcoLine F/Optic

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

• \$ 979.00 (F/Optic)

• \$ 693.00 (N/Optic)



1 Year FairCare Warranty
Unconditional warranty except physical damage

Buy 3 - Get 10% OFF

ECO LINE POWER EDITION



- 28 Watt Power
- 60 dB Operating Noise
- AeroPerformance blade
- Glass rod fiber optics
- Refined ceramic bearings
- Smooth Eco body

• HE20KL
\$1,265.00

Buy 3 or more
Only \$1,089 ea



With 18 months
FairCare Warranty
Unconditional warranty
except physical damage

Available for only
KaVo Multiflex coupling fitting

ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 - Standard head with 4 hole water spray
- HE18 - Small head with 3 hole water spray



\$ 539.00

Buy 3 - Get 10% OFF

**Buy 3 - Get 20% OFF
Mix & Match**

**NEW
RANGE**

**With 2 Years
FairCare Warranty**





**SAFE-TECH
PUSH
BUTTON**



 1:5 Red ●


 Max. 200,000 rpm

 Internal water


 For FG-burs
(Ø 1.6 mm)


● LP15L (F/Optic)
RRP: \$1,826.00

Buy 3 - \$1,460.80 ea

 1:1 Blue ●

 Max. 40,000 rpm


 Internal water

 For CA-burs
(Ø 2.35 mm)


● LP11L (F/Optic)
RRP: \$1,287.00

Buy 3 - \$1,029.60 ea

 5.4:1 Green ●

 Max. 7,407 rpm

 Internal water


 For CA-burs
(Ø 2.35 mm)


● LP54L (F/Optic)
RRP: \$1,320.00

Buy 3 - \$1,056.00 ea

 1:1 Blue ●

 Max. 40,000 rpm

 Internal water

 For HP-burs
(Ø 2.35 mm)

● LP01L (F/Optic)
RRP: \$1,045.00

Buy 3 - \$836.00 ea

AIR MOTORS 



NEW
Shorter
and
Lighter

- ✓ Internal water
- ✓ LED

● AM1006
\$ 990.00



- ✓ Internal water

● AM1004
\$ 748.00



- ✓ External water

● AM0014
\$ 473.00



● AM0004
\$ 473.00

NEW

KAVO COMPATIBLE HEADS 



- Kavo 68LH compatible
- Blue (1:1)
- LCH11
\$ 473.00



- Kavo 67LH compatible
- Green (2:1)
- LCH21
\$ 473.00



- Kavo 31LR compatible
- Green (2:1)
- LCH21P
\$ 484.00

Buy 3 - Get 1 FREE

- Mix & Match
- Cheapest Free



FairCare Warranty¹
1 year unconditional warranty
except physical damage



⚙️ 1:5 Red ●
○ Internal water
For FG-burs
(Ø 1.6 mm)

• LE15L (F/Optic)
RRP: \$1,595.00

Buy 3 - Get 1 FREE

⚙️ 1:1 Blue ●
○ Internal water
For CA-burs
(Ø 2.35 mm)

• LE11L (F/Optic)
RRP: \$1,155.00

Buy 3 - Get 1 FREE

⚙️ 5.4:1 Green ●
○ Internal water
Screw-in & Snap-on
Prophy Cups

• LE54LP (F/Optic)
RRP: \$1,287.00

Buy 3 - Get 1 FREE

⚙️ 1:1 Blue ●
○ Internal water
For HP-burs
(Ø 2.35 mm)

• LE01L (F/Optic)
RRP: \$990.00

Buy 3 - Get 1 FREE

• LE15 (Non-optic)
RRP: \$1,375.00

Buy 3 - Get 1 FREE

• LE11 (Non-optic)
RRP: \$891.00

Buy 3 - Get 1 FREE

• LE54P (Non-optic)
RRP: \$1,045.00

Buy 3 - Get 1 FREE

• LE01 (Non-optic)
RRP: \$825.00

Buy 3 - Get 1 FREE

Mix & Match



**Buy More
Save More**



⚙️ 1:1 Blue ●
○ Push button
• LB11/11
Buy 1 - \$418.00

Buy 3: \$339.00 ea

Buy 4+ \$319.00 ea

⚙️ 1:1 Blue ●
○ Latch Type
• LB11/11C
Buy 1 - \$418.00

Buy 3: \$339.00 ea

Buy 4+ \$319.00 ea

⚙️ 4:1 Green ●
○ Snap-on Cups
• LB41/11P
Buy 1 - \$594.00

Buy 3: \$484.00 ea

Buy 4+ \$473.00 ea

⚙️ 1:1 Blue ●
○ For HP-burs
• LB02
Buy 1 - \$418.00

Buy 3: \$339.00 ea

Buy 4+ \$319.00 ea

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019



EMS Non-optic
Type EN-041

- SC21E
Buy 1 - \$462.00

Buy 3 - Get 1 FREE



EMS Piezon Non-optic
Type EN-061

- SC21E2
Buy 1 - \$561.00

Buy 3 - Get 1 FREE



EMS Piezon LED
Type EN-060A

- SC21EL
Buy 1 - \$748.00

Buy 3 - Get 1 FREE



KaVo Piezo LED
Type 1.007.3995

- SC21KL* **NEW**
Buy 1 - \$748.00

Buy 3 - Get 1 FREE



Acteon-Satelec Non-optic
Type F 12 281

- SC21S
Buy 1 - \$462.00

Buy 3 - Get 1 FREE



Acteon-Satelec LED
Type F 12 609

- SC21SL
Buy 1 - \$748.00

Buy 3 - Get 1 FREE

*Only for use with EMS compatible scaler tips


▼ SCALER TIPS




Buy 3 - Get 1 FREE
while stock lasts

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface


• **SCALING DLC COATED** Only available for Satelec fitting

 ST03SD
General Scaling
\$88.00

 ST04SD
SATALEC 2
\$132.00

• **PERIO DLC COATED** Only available for Satelec fitting


 ST06SD
SATELEC 10Z
\$132.00

 ST10SD
SATALEC 10P
\$154.00

• **ENDO DLC**


 ST08SD
SATALEC S07
\$154.00


• **SCALING STAINLESS STEEL** Only available for Satelec fitting

 ST03S
General Scaling
\$63.80


 ST04S
SATELEC 2
\$88.00

• **PERIO STAINLESS STEEL** Available for EMS & Satelec fittings

 ST07E / ST07S
EMS DPL3
SATELEC 10X
\$110.00

 ST09E / ST09S
EMS HPL
\$110.00

 ST10E / ST10S
SATALEC 10P
\$110.00

 ST11E / ST11S
EMS RT1
\$110.00

Universal Torque Wrench



TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening




Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



\$1,210.00


Lightweight
ONLY 135 gms

Order code: 15501

3 Easy Steps Process



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



www.anthogyr.com

Anthogyr
A Straumann Group Brand

▼ PROPHY LINE AIR POLISHING DEVICE

**Buy 1 Polisher and receive
• 1 bottle of Prophy Powder FREE
and
• 50% discount on 1 additional tip**



- **PR1011**
\$1,375.00
- Ergonomically shaped, compact and light-weight air-driven polishing device

- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation

- **Prophy Powder (300gms) Strawberry - PP1013**
\$ 60.50 each

Additional tip



- **SP3011**
\$ 329.00

Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

▼ LEDs

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



\$88.00

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

▼ PREMIUM LUBRICATION OIL



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors
- **LU1011 (500 ml)**
\$ 68.20



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system
- **LU1022 (500 ml)**
\$ 68.20



- ✓ For W&H Assistina lubrication system
- **LU1031 (500 ml)**
\$ 88.00

MANUFACTURER



MK-dent GmbH
Marie-Curie-Str. 2
22941 Bargteheide
GERMANY
www.mk-dent.com

YOUR SPECIALIST DEALER



• All catalogue prices are inclusive of GST • Promotion valid until 31st Dec 2024 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.